

Wearanize+ Participant To-Do List

1. Please fill out the **Evening questionnaires** using the link received via email.

Before going to bed:

2. Turn on the **Zmax headband** by flicking the SWITCH (1) and pressing the MIDDLE BLACK BUTTON (2). It will make sounds as you press them. If the recording has started, the blue lights under the device will turn off.
3. Turn on the **Somnoscreen** by pressing the ON (1) button once and the ENTER (4) button **four** times.

Actions: 1: Turn on > wait 5 seconds >

4: Record > 4: Ok > 4: Ok > 4: Start > Recording.

If successful, it will show “Recording...” and the battery percentage. If not, please follow the menu shown in the monitor and start over.

See a video using the QR code for more help:



4. To set up the **Gmeter**, attach the headband to the Gmeter box (do not use the screws).

5. Then, connect the accelerometer to the **Gmeter** box.

6. To turn on the **Gmeter**, flick the “ON/OFF” switch. See the **red square** in the picture. Then flick the “REC” switch. See the **green square** in the picture.

If the recording has started, all three LED lights on the device should turn blue. If it did not work, please repeat the steps again.



7. Please jump **five** times, as practiced in the lab. Please remember to hold all the devices/boxes in place with your hands while jumping.
8. Go to bed.
9. While in bed, perform these activities before falling asleep (in this order):

- a. close your eyes and relax for 1 minute
- b. with your eyes closed, clench your teeth for 2 seconds, then let go for 2 seconds. Repeat 5 times.
- c. with your eyes closed, perform the Left-Right (LR) eye movement 5 times.
- d. blink 5 times.

Checklist:

- ☐ Jump: 5 times
- ☐ Close eyes & relax: 1 min
- ☐ Teeth clench: 5 times
- ☐ LRLR eye movement: 5 times
- ☐ Blink: 5 times

10. Try to fall asleep.

After waking up:

1. After waking up, perform these activities while staying in bed (in this order):
 - a. close your eyes and relax for 1 minute
 - b. with your eyes closed, clench your teeth for 2 seconds, then let go for 2 seconds. Repeat 5 times.
 - c. with your eyes closed, perform the Left-Right (LR) eye movement 5 times.
 - d. blink 5 times.

Checklist:

- ☐ Close eyes & relax: 1 min
- ☐ Teeth clench: 5 times
- ☐ LRLR eye movement: 5 times
- ☐ Blink: 5 times

2. Leave the bed.
3. Jump **five** times again. Please remember to hold all the devices/boxes in place with your hands while jumping.

4. Please take off **both wristwatches**. Turn off the **Empatica** (the wristband with no display) by pressing the button for one second.
5. Take off the **Zmax headband** by loosening the straps and turn it off by flicking the SWITCH (1) (see previous page)
6. To turn off the **Gmeter**, first, flick the “REC” switch (see the **green square** on the previous page), then flick the “ON/OFF” switch (see the **red square**). Please flick the switches in the mentioned order. Otherwise, the recording will not be saved.
7. Take off the **Gmeter headband** from your forehead. Do not remove the tape from the headband.
8. Detach the **Gmeter headband** and the **accelerometer** from the **Gmeter box**.
9. Take off the cap and the (wired) electrodes one by one. Please do not pull too hard. Check if the contact with the skin is loose if they do not come off easily.
10. Take off the **Somnoscreen** from your waist.
11. Turn off the Somnoscreen by pressing three buttons (2, 3, 4) together. Then press 3, 4, 3, 4. The display will turn off if successful.
 Actions: 234:Stop recording > 3: to Cancel > 4:Cancel > 3:to Sure > 4:Sure.
12. Put the gluey ends of all the electrodes in the zip lock bag (provided). Please do not detach the electrodes from the box.
13. Take off the **ActivPAL** leg patch.
14. Wash your hair with shampoo (a bit more than usual and for a bit longer).
15. Fill out the **Morning questionnaire** using the link received via email.
16. Put the **Gmeter** in the **small box** and the other devices in the **big box/the bag**.
17. Drop off the bag with all the devices/boxes to the discussed place before noon.



Somnoscreen:



Zmax:



mSafety:



Empatica:



ActivPAL:



Gmeter:



For any questions/problems with the devices/recordings, please contact via email: niloy.sikder@donders.ru.nl or phone/text: +4917679021503