

After a participant signs up through Limesurvey/Sona:

1. Contact them with informed consent form, other general forms, and the procedures to open a Sona account (if necessary).
2. Fix a date and time for the session
Preliminary: wear comfortable clothing, no jewelry, tie your hair in a (high) ponytail (if long hair), have an early dinner (if possible) or have a light dinner afterwards, try to go to sleep soon after the session.
3. Find a buddy and/or a lab member with BROK
4. Book the lab
5. Send reminders (1 week before, the day before)
Once the participant is in the lab:
6. Check their suitability (head size and sleeping position)

If not suitable: They get paid a certain amount if they sign the forms

If suitable:

7. Explain the info on the informed consent and other forms (regulations and rights). They can abort at any time without reason.
8. Explain the study procedure
9. Get the informed consent form signed
10. Screening in Castor.
11. Explain what they have to do after the session (participant checklist).
12. Decide on a device pickup location, time, and person.
13. They fill out the screening form in Castor.
14. Initialize the SD cards of Somnoscreen and Zmax and check the charge parentages of the other devices.
15. Change the Somnoscreen battery (with a fully recharged one)
16. Put the devices/electrodes on in this order:
 - a. ActivPAL(1): show them the position, put some tape on it, they go to the washroom/another room and put it on themselves
 - b. Empatica(2) (non-dominant hand) and Msafety(3) (dominant hand)
 - c. Somnoscreen(4) (ECG electrodes(4a): right shoulder and left belly), ask if they want to do it themselves
 - d. Somnoscreen (position the box(4b) on the chest)
 - e. Check hair
 - f. Measure the head with tape to determine cap size, Cz position
 - g. Put the cap on and mark: C3, C4, F3, F4, O1, O2, Gnd (check if the touches are too hard)
 - h. Mark/identify the Mastoid and EOG positions
 - i. Put the EEGs(4c) on, one by one (check if the touches are too hard)

- j. Before putting the EOGs on, put the Zmax (and Gmeter combo) on the forehead and check how much space it takes. If they overlap, adjust the position.
 - k. Put the EOGs(4d) on.
 - l. Put the EMGs(4e) on (use tape).
17. Show them how to turn the Somnoscreen recording on (and off).
 18. Clean the forehead with a patch.
 19. Show them the Zmax buttons and how to turn on, record, and turn off
 20. Put the Zmax(5) and Gmeter(6) combo on, don't turn on the recording.
 21. Put on a net/cap on the head, and check for any loose connections
 22. Tell them about the synchronization movements (jumps and LR eye movements)
 23. Show them how to jump
 24. Tell them to stand up and jump 5 times. Check the height and frequency.
 25. Check for any loose connections again. If so, fix them
 26. Explain to them how to take the devices off the body in the morning. Show them the video.
 27. Hand them a bag/box to store the devices, a zip lock bag for the electrodes
 28. Hand them the participant checklist, tell them to read it, and ask if they have any questions
 29. Send them the questionnaires from Castor, 3 immediately, 2 scheduled in the morning
 30. Tell them to fill in the questionnaires: three at night, two in the morning, separate links.
 31. Remind them of the device pickup location, time, and person.
 32. Put the Somnoscreen battery in charge.
 33. Lab logbook entry
 34. Finish screening in Castor
 35. Send the questionnaires.
 36. Clean the cap and clean up the lab.

After receiving the devices:

1. Extract/copy data to P:\\3028005.01 drive using softwares:
 - a. Somnoscreen: Dominos
 - b. Zmax: HDrecorder
 - c. Empatica: E4 Manager
 - d. Activpal: PALconnect
 - e. Gmeter: WearableGlymphometer
2. Recharge all devices
3. Clean electrodes:
 - a. wear gloves,
 - b. prepare solution (50 mL Sekusept + 950 mL water),
 - c. clean the contact ends of the electrodes using a brush with hot water in a bowl/bucket
 - d. Soak the contact ends (and the Cap) into the solution for 15 mins (set an alarm).
 - e. Rinse with tap/distilled water
 - f. Hang the cap to dry
 - g. Let the electrodes dry on a towel